


# Wall Mounting Cabinets

LabFit's Wall Mounting Cabinets feature adjustable shelving and optional glass viewing windows. All wall cabinets utilize a French cleat system for a low-profile and simple to install wall mount.

## With Glass Windows

Part #	Description	Dimensions WHD
53699	Cabinet Wall Single Width White Polypropylene One Left Door with Glass	24" X 30" X 13"
53698	Cabinet Wall Single Width White Polypropylene One Right Door with Glass	24" X 30" X 13"
53700	Cabinet Wall Double Width White Polypropylene Double Door with Glass	30" X 30" X 13"
53701	Cabinet Wall Double Width White Polypropylene Double Door with Glass	36" X 30" X 13"
53702	Cabinet Wall Double Width White Polypropylene Double Door with Glass	46" X 30" X 13"



For pricing and quotes on LabFit products, contact the TrippNT Sales team:

Email: [Sales@TrippNT.com](mailto:Sales@TrippNT.com)

Phone: 816-792-2604



Industry-Leading  
Warranty



Fast Lead Times

