


Wall Mounting Cabinets

LabFit's Wall Mounting Cabinets feature adjustable shelving and optional glass viewing windows. All wall cabinets utilize a French cleat system for a low-profile and simple to install wall mount.

Standard Doors

Part #	Description	Dimensions WHD
53656	Cabinet Wall Single Width White Polypropylene Right Hand Door	18" X 30" X 13"
53657	Cabinet Wall Single Width White Polypropylene Left Hand Door	18" X 30" X 13"
53662	Cabinet Wall Single Width White Polypropylene Right Hand Door	24" X 30" X 13"
53703	Cabinet Wall Double Width White Polypropylene Double Door	30" X 30" X 13"
53704	Cabinet Wall Double Width White Polypropylene Double Door	36" X 30" X 13"
53661	Cabinet Wall Double Width White Polypropylene Two Doors	48" X 30" X 13"



For pricing and quotes on LabFit products, contact the TrippNT Sales team:

Email: Sales@TrippNT.com

Phone: 816-792-2604



Corrosion Resistant



Easy to Install

