



LabFlow Trays

Polypropylene



Polycarbonate



MODULAR TRAYS THAT MOVE WITH YOUR WORKFLOW

LabFlow Trays are the flexible foundation of the entire LabFlow storage system. Designed for compatibility across all LabFlow carts and racks, these high-visibility trays simplify organization, speed up supply access, and let you re-configure your workspace instantly. Whether you're storing reagents, experiments, or daily-use tools, LabFlow Trays keep your lab running clean, contained, and in control.

MODULAR BY DESIGN

- Compatible with all LabFlow carts and racks for seamless, cross-system organization.
- Easily swap, reorganize, or transport trays to fit changing workflows.

THREE DEPTH OPTIONS


- Available in 2.5", 4", and 6" depths to optimize vertical space and Storage needs.

HIGH VISIBILITY

- Translucent or clear construction makes contents instantly identifiable—no opening or guessing.

LIFETIME GUARANTEE


- Work stress-free, knowing your trays are covered for life.



**Compatible
with any
LabFlow
Cart!**



LabFlow Trays



Thirteen 4" Trays

Nine 6" Trays

Six 4" Trays
Five 6" Trays

Which Material to Choose?

Polypropylene

- Translucent color
- Cheaper and lighter weight
- -40°F to 160°F temperature range
- BPA free and chemical-free
- NSF listed

Polycarbonate

- Totally clear
- Virtually unbreakable
- -40°F to 212°F temperature range
- NSF listed



Corrosion Resistant



Short Lead Times




Ships Assembled




Made in the U.S.A



100% Lifetime Guarantee




Three 6" Trays




Three 4" Trays



Six 2.5" Trays




Two 2.5" Trays
Two 6" Trays



Three 2.5" Trays
Two 4" Trays

LabFlow Trays



		Internal Dimensions		
Material	Part Number	Width	Length	Depth
Polypropylene	101171	10 ¾"	18 ¾"	2 ½"
	101126	10 ½"	18 ½"	3 ¾"
	101170	10 ¼"	18 ¼"	5 ¾"
Polycarbonate	101172	9 ½"	17 ¼"	2 ¼"
	101167	10"	17 ¾"	3 ¾"
	101120	10"	18"	5 ¾"