



Double ABS Bottle Holder 500mL

Part No. 50017


TrippNT's Double ABS Bottle Holder 500mL keeps your bottles upright and secure. 3 1/8 inch Diameter Holes 7 1/2 X 3 1/4 X 4 1/8 inches WHD

Rev: 4



Features:

- Open Rectangular Design with Top Circle Slot
- Mounts on Wall with Pre-Drilled Holes
- Constructed from White ABS
- Holds 2 500mL Bottles



Made in the U.S.A



Corrosion Resistant



100% Lifetime Guarantee