






Shallow 4 Slot Drawer Organizer

Part No. 50040

Keep your most-used lab tools organized and within easy reach with the TrippNT Shallow 4-Slot Drawer Organizer Tray. Designed for lab managers who value efficiency and order, this compact tray fits neatly into drawers or on benchtops to separate and protect delicate instruments, pipettes, and small supplies. Made from durable white styrene, it's lightweight, easy to clean, and built to withstand the demands of daily lab use—helping your team work faster and stay focused on what matters most.

13 5/8" x 1 3/4" x 8 7/8" inches WHD



Features:

- Ideal for pipettes, markers, spatulas, scissors, or other everyday items.
- Use inside drawers or on benchtops for flexible organization.
- Smooth, nonporous surfaces clean easily with common lab disinfectants.



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee