



Large Draft Shield





Part No. 50046

When a gust of air distorts crucial measurements, bad data costs time and trust. This shield isolates your instrument, preserving accuracy and confidence so your experiments don't have to be second-guessable.

18 X 18 1/8 X 19 1/4 inches WHD

Features:

- Medium Cube-Designed Shield with Draft Deflecting Overhang
- Semi Circle Cut Out for Cord and Wire Passthrough
- Constructed from Clear 3/16" Acrylic
- Blocks Air Drafts around Delicate Equipment



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee

