



16 Compartment Drawer Organizer


Part No. 50060

In a multi-lab facility seconds count. The 16-Compartment Drawer Organizer offers fourteen compartments so you can sort vials, tips, and small consumables with clarity and speed.

17 1/2 x 2 3/8 x 17 1/2 inches WHD

Features:

- 16 small compartments.
- Rounded inside corners allow easy cleaning.
- Great for separating small consumables or tools.
- Can be used as a tray or drawer organizer



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee