



14 Compartment Drawer Organizer

Part No. 50063


In a multi-lab facility seconds count. The 14-Compartment Drawer Organizer offers fourteen compartments so you can sort vials, tips, and small consumables with clarity and speed.

19 x 2 $\frac{3}{8}$ x 17 $\frac{3}{8}$ inches WHD



Features:

- 14 small compartments.
- Rounded inside corners allow easy cleaning.
- Great for spill containment.
- Can be used as a tray or drawer organizer



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee