






## Lab Storage Shelf with 4 Bins

Part No. 50065

Lab techs want clutter-free benches—but small supplies often get scattered, wasting time and focus. The Lab Storage Shelf keeps essentials organized and visible, leading to faster access and tidier workflows.

24 x 12 ½ x 10 ⅛ inches WHD



### Features:

- 4 Storage Bins with Top Shelf
- 1.5-inch Safety Lip
- Constructed from White PVC
- Stores Frequently Used Supplies in Divided Bins



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**