





## 4 Compartment Drawer Organizer

Part No. 50066

With four equal compartments, this white styrene organizer sorts vials, pipette tips, labels, and tools cleanly. Use it as a tray or drawer insert.

16 <sup>3</sup>/<sub>4</sub> x 2 <sup>3</sup>/<sub>8</sub> x 12 <sup>3</sup>/<sub>4</sub> inches WHD



### Features:

- 4 compartments.
- Rounded inside corners allow easy cleaning
- Great for small consumables and tools
- Can be used as a tray or drawer organizer



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**