



## Large Pipette Storage Tray with Adjustable Dividers

Part No. 50069


Store, separate, and protect your pipettes the smart way. This spacious PVC tray features 4 adjustable dividers to accommodate a wide range of pipette sizes and accessories. Designed to fit neatly inside drawers or sit on shelves, it's a durable, chemical-resistant solution for labs that need reliable organization without wasted space.

25 1/8 x 3 1/4 x 15 3/4 inches WHD



### Features:

- **Flexible Storage:** 4 movable dividers let you customize spacing for different pipettes and tools
- **Drawer- or Shelf-Ready:** Low-profile design fits most standard lab drawers
- **Protects Instruments:** Prevents damage and clutter by keeping pipettes separated and secure
- **Durable PVC Build:** Withstands daily lab use and harsh cleaning agents



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**

