



## Medium Draft Shield





Part No. 50073

When a gust of air distorts crucial measurements, bad data costs time and trust. This shield isolates your instrument, preserving accuracy and confidence so your experiments don't have to be second-guessable.

15 X 15 X 16 1/4 inches WHD

### Features:

- Medium Cube-Designed Shield with Draft Deflecting Overhang
- Semi Circle Cut Out for Cord and Wire Passthrough
- Constructed from Clear 3/16" Acrylic
- Blocks Air Drafts around Delicate Equipment



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**