



## Single ABS Bottle Holder 1000 mL


Part No. 50080

In busy labs, cluttered benches and unstable solvent bottles can create hazards and disrupt workflow. This wall-mounted ABS bottle holder gives you back your space while keeping 1 L bottles upright and accessible. Its open design accommodates low-volume bottles and reduces fume exposure from tipped containers.

4  $\frac{3}{4}$  x 4  $\frac{1}{4}$  x 4  $\frac{7}{8}$  inches WHD

### Features:

- Holds one 1000 mL bottles in 3  $\frac{3}{4}$  diameter circle slots
- Mounts directly to the wall to free up crowded benchtop space
- Open top and front provide easy access to bottles at all times
- Durable and chemical-resistant easy to clean with standard disinfectants



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee