



## Safety Shelf Workstation

Part No. 50104




Organize your lab bench with this extensive workstation. Keep clutter out of the way and create more space to do more meaningful work. The safety lip acts as protection and a barrier for your benchtop against leaks and spills.

22 7/8 x 14 x 9 7/8 inches WHD



### Features:

- Consists of 3 shelves, 2 side compartments, and 2 tool holders.
- 1 Inch safety lip on the three shelves to protect from spills or leaks.
- Constructed from easy to clean PVC



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**