




## Straight Triple Safety Shelf

Part No. 50113

Protect your supplies with TrippNT's Triple Safety Shelf. The workstation's lipped design protects your bench and surrounding area from any leaks or possible stains while keeping your benchtop organized.

23 7/8 x 6 5/8 x 16 inches WHD



### Features:

- 3 Shelves
- 1 1/2 inch Safety Lip
- Protects Your Bench and maintains organization
- Constructed from white PVC



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee