



## Large Step Shelf 24 inch

Part No. 50123


The perfect shelf to use right at your work area to add space. This shelf does not attach to a wall, but stands freely. This allows it to be moved around exactly where you need it, when you need it.

23 <sup>3</sup>/<sub>4</sub> x 10 <sup>1</sup>/<sub>8</sub> x 8 <sup>5</sup>/<sub>8</sub> inches WHD



### Features:

- Raised 10in to fit over other bench ware
- Cut out across the bottom for cord pass through
- Constructed with PVC making this shelf easy to clean



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**