



## Lab Storage Shelf 2 shelf & 9 Bins


Part No. 50152

Organize critical supplies with confidence. TrippNT's Medium Lab Storage Shelf creates flexible, adjustable storage for bottles, gloves, wipers, or chemicals, all in one durable, easy-to-clean unit. With two sturdy shelves and seven removable dividers.

12 x 20 1/8 x 10 1/4 inches WHD

### Features:

- Two shelves and seven removable dividers create nine customizable bins for bottles, gloves, or reagents.
- Made from cleanable, chemical-resistant white PVC
- Weighs under 5 lbs
- Work stress-free, knowing your product is covered for life.



Unit ships with 3 large dividers, and 4 small dividers.



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**