



Quadruple ABS Bottle Holder 500mL

Part No. 50157


Lab and clinical staff need safe, clean bottle storage—without spills or clutter. TrippNT's Quadruple ABS Bottle Holder secures four 500 mL bottles in an upright, accessible position for a tidier, more efficient workspace.

14 $\frac{7}{8}$ x 3 $\frac{1}{4}$ x 4 $\frac{1}{8}$ inches WHD



Features:

- Open Rectangular Design with Top Circle Slot
- Mounts on Wall with Pre-Drilled Holes
- Constructed from White ABS
- Holds Four 500mL Bottles



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee