






Fume Hood Shelf 3 Feet

Part No. 50215

Designed specifically for fume hoods and contamination boxes, TrippNT's 3 Foot Fume Hood Shelf offers a stable, contamination-conscious surface for essential supplies. It's made from white PVC, a cleanable, chemical-resistant material, and allows unrestricted laminar flow to protect the integrity of your hood.

36 x 8 x 5 7/8 inches WHD



Features:

- Designed to maintain proper air pull inside the hood for safe operations
- Easy to Clean
Chemical-Resistant PVC
- Keeps critical supplies within reach while preserving contamination control
- 1-inch edge keeps items securely in place, preventing spills or drops inside the hood



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee