



## HDPE Fume Hood Shelf 3 Feet




Part No. 50315

Designed specifically for fume hoods and contamination boxes, TrippNT's 3 Foot Fume Hood Shelf offers a stable, contamination conscious surface for essential supplies. It's made from white HDPE, a cleanable, chemical-resistant material, and allows unrestricted laminar flow to protect the integrity of your hood.

36 x 8 x 5 7/8 inches WHD

### Features:

- Designed to maintain proper air pull inside the hood for safe operations
- Easy to Clean  
Chemical-Resistant HDPE
- Keeps critical supplies within reach while preserving contamination control
- 1-inch edge keeps items securely in place, preventing spills or drops inside the hood



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**