



## Extra Large Bench Top Workstation

Part No. 50377





Keep your workstation free of clutter with TrippNT's Extra Large Benchtop Organizer. Keep consumables and other supplies exactly where you need them, when you need them.

10 x 30 1/8 x 15 7/8 inches WHD



### Features:

- Store Bottles, Gloves, Wipers, or Chemicals right at your fingertips
- Made with 4 shelves with 15 total compartments and 6 tool holder slots on the back
- Constructed with white PVC which makes this easy to clean



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**