



Quadruple HDPE Bottle Holder 500 mL

Part No. 50539


In busy labs, cluttered benches and unstable solvent bottles can create hazards and disrupt workflow. This wall-mounted HDPE bottle holder gives you back your space while keeping 500 mL bottles upright and accessible. Its open design accommodates low-volume bottles and reduces fume exposure from tipped containers.

14 ⁵/₈ x 3 ¹/₈ x 4 ¹/₈ inches WHD



Features:

- Holds four 500 mL bottles in 3 ¹/₈ inch diameter circle slots
- Mounts directly to the wall to free up crowded benchtop space
- Open top and front provide easy access to bottles at all times
- Durable and chemical-resistant HDPE easy to clean with standard disinfectants



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee

