



## Triple HDPE Bottle Holder 500 mL

Part No. 50571


In busy labs, cluttered benches and unstable solvent bottles can create hazards and disrupt workflow. This wall-mounted HDPE bottle holder gives you back your space while keeping 500 mL bottles upright and accessible. Its open design accommodates low-volume bottles and reduces fume exposure from tipped containers.

11 1/8" x 3 1/8" x 4 1/8" inches WHD



### Features:

- Holds three 500 mL bottles in 3 1/8 inch diameter circle slots
- Mounts directly to the wall to free up crowded benchtop space
- Open top and front provide easy access to bottles at all times
- Durable and chemical-resistant HDPE easy to clean with standard disinfectants



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**

