



# Lab Islands Series



## MAXIMUM MOBILE WORKSPACE. TOTAL LAB FLEXIBILITY.

Every lab professional wants a workspace that adapts as quickly as their needs change. But most lab furniture is heavy, fixed in place, or hard to clean—making it a challenge to reconfigure spaces, secure supplies, or maintain a sterile environment. That’s frustrating when efficiency and flexibility are non-negotiable. TrippNT’s Lab Island Carts solve this by combining cleanable HDPE construction, secure lockable storage, and smooth mobility in one versatile unit. You get a reliable, movable workstation that makes your lab cleaner, safer, and easier to work in—day after day.

### EASY MOBILITY

- Heavy-duty casters make it simple to move or rearrange units within labs, classrooms, or between spaces.

### ROBUST WORK SURFACE

- HDPE top gives you a span of clean, resistant workspace for prep, equipment, or instruments.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

## ASK ABOUT OUR DEMO PROGRAM

Try it out at your facility.



#### Full Extension Shelf

- Retrieve supplies easily without any strain.



#### Easy and Secure Mobility

- Casters make maneuvering the Lab Island easy; lock them to secure it in place.



#### Lockable Tambour Door

- Safely store all supplies or equipment when needed.



#### Large & Durable Worktop

- The HDPE worktop is durable, stain- and corrosion-resistant.



## COLOR OPTIONS

Belize



PN: 50593  
White

Cayman



PN: 50592  
White

Maui



PN: 50591  
White



PN: 52943  
Black



PN: 52937  
Black



PN: 51002  
Black

Model	Width	Height	Depth	Weight
Maui	24 ¼"	35 ⅜"	23 ⅞"	79 lb
Cayman	47 ¾"	35 ½"	23 ⅞"	180 lb
Belize	93 ½"	35 ½"	23 ⅞"	275 lb



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee