





28 Compartment Adjustable Reagent Tray

Part No. 50712

Maximize lab reagent access and transport with this adjustable tray that supports up to 28 bottles. Built from white PVC, with four built-in handles for safe movement around busy cold-rooms and prep areas

18 7/8 x 11 1/4 x 2 3/4 inches WHD



Features:

- 28 compartments.
- Removable dividers allow easy cleaning
- Great for transporting reagents or samples
- Can be used as a tray or drawer organizer



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee