



# MRI Narrow Lockable Cart With Roll Top Door

50729



## STRENGTH BEHIND THE SCENES

In the MRI suite, every detail counts, and every moment demands focus. The TrippNT MRI Narrow Locking Cart doesn't compete for your attention, it supports it. Its slim frame and tambour door allows you to tuck it away securely and move it freely right when you need it. Built with MRI-safe materials, it blends seamlessly into the environment, protecting both staff and patients. With the quiet assurance that supplies are secure, accessible, and exactly where they're needed, caregivers can stay centered on the work that matters most.

## MRI-SAFE BUILD

- Non-ferrous materials and casters ensure compatibility and safety in magnetic environments.

## ROLL TOP SECURITY

- Sliding tambour door locks supplies away without blocking access or hallways

## LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

## [ASK ABOUT OUR DEMO PROGRAM](#)

*Try it out at your facility.*



### Ergonomic Handle

- Transport supplies with ease and control



### Removable Tray System

- Restock and customize with ease



### Tambour Door

- Secures and protects supplies without taking up extra space.



### Locking Hasp

- Durable plastic, lock-compatible side hasp for added security.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | [Sales@TrippNT.com](mailto:Sales@TrippNT.com)

# MRI Narrow Lockable Cart With Roll Top Door

50729

## EASY DISINFECTION

- Non-porous HDPE and Styrene surfaces clean quickly with standard hospital disinfectants.

## CUSTOMIZABLE STORAGE

- Six pocket trays adjust to changing supply needs, supporting flexible workflows.

## CONTROLLED MOBILITY

- Combination of swivel and locking casters keeps the cart steady when parked and smooth in motion.

## QUIET OPERATION

- Smooth door and tray movement minimize noise, preserving a calm environment for patients.

## SLIM DESIGN

- Fits in tight MRI suites while still offering ample storage capacity.



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**







50424 Tray Dimensions:

8 3/8" Width

2 3/4" Height

12 1/2" Depth



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com