



MRI Belize Lab Island

50734



Maximum Workspace & Coil Storage for MRI Suites

In a busy MRI suite filled with tight corridors, heavy coils, and constant patient turnover, MRI technologists need equipment that works *with* them — not against them. The TrippNT MRI Belize Lab Island Cart delivers exactly that: an 8-foot, MRI-safe mobile workstation with lockable roll-top storage and wide shelving built for speed, safety, and organization. Its non-ferrous HDPE construction, smooth-rolling casters, and generous workspace make it easy to stage coils, swap linens, and access supplies without stepping away from the scan room. Designed for technologists who thrive under pressure, it keeps workflows efficient, calm, and compliant.

100% MRI-CONDITIONAL MATERIALS

- 99.4% Non-Metal construction. Zero ferrous metals.

DURABLE, CHEMICAL-RESISTANT SURFACES

- HDPE construction resists stains, disinfectants, and daily wear, preserving aesthetics.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



[ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



LOCKING MECHANISM FOR SECURITY

- Integrated lock ensures spare linens, coils, or sensitive supplies stay in place when unattended.



INTERNAL STATIONARY SHELF

- Keeps frequently used items within reach without needing to dig.



2 SLIDING TAMBOR DOORS

- Secures gear with a roll-up door that disappears when open—no swinging panels in tight space.



SMOOTH, LARGE CASTER

- 5" twin-wheel casters (3 locking, 3 swivel) provide easy gliding over scan room floors.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

MRI Belize Lab Island

50734

HIGH CAPACITY MODULAR STORAGE

- shelf supports coil placement during prep or change-over, keeping equipment stable while freeing both hands for patient setup.

KEEP COILS AND PATIENT POSITIONERS CLOSE

- Stay on schedule by creating a flexible and mobile workflows

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Corrosion Resistant



Short Lead Times






Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



RELATED PRODUCTS



MRI Maui

51820



MRI CAYMAN

51821