



Procedure Cart Series



Engineered for Urgency, Built for Reliability

In fast-paced clinical settings, the TrippNT 51046 Procedure Cart keeps anesthesia teams prepared with secure, clutter-free readiness. Six locking drawers, a recessed tray top, and smooth-gliding casters streamline organization and mobility. Built in the USA from non-metal HDPE, it arrives fully assembled—delivering dependable performance in tight, high-stakes clinical spaces.

Hygienic & Durable

- Corrosion-resistant HDPE and ACP construction resists disinfectants and avoids rust, safeguarding equipment integrity

Compact Yet Capacity-Rich

- Wide-push design packs space-saving storage for anesthesia kits without crowding tight OR corridors

Secure, Mobile Readiness

- Six full-extension, key-locking drawers and 4-inch twin-wheel casters (2 locking, 2 swivel) keep instruments both organized and locked during transport.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

ASK ABOUT OUR DEMO PROGRAM

Try it out at your facility.



Recessed Top

- Provides a stable, secure work surface



Wide Drawers

- Fits up to two 11in trays in each drawer



Two Drawer Depths

- Large capacity for a range of equipment



Narrow Footprint

- 20" Depth enables use in tight spaces



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

Procedure Cart Series

Streamlined Workflow




- Built-in recessed tray top with easy-grab handle gives clinicians quick-access staging for tools during procedures.

Ready Out of the Box

- Ships fully assembled, eliminating setup delays and ensuring consistent, reliable quality from day one.

Smooth Operation

- Balanced 88 lb. design and 4-inch casters glide quietly and steadily through crowded corridors without jostling supplies.



COLOR OPTIONS



PN: 51046
White



PN: 51078
Almond Beige



PN: 51071
Global Blue



PN: 51358
Black



PN: 53593
Dolphin Gray



PN: 51079
Silver Metallic

ACCESSORIES



PN: 50917
4 Compartment



PN: 50974
6 Compartment



PN: 50975
4 Compartment



PN: 51998
Adjustable