



Double Two Face Glove Box



Part No. 51144

Enhance organization and accessibility in your shared lab spaces with the TrippNT Two-Faced Double Glove Box Holder. Designed for lab managers who value efficiency and cleanliness, this dual-access dispenser allows glove retrieval from both sides of a pass-through wall or shared workstation—keeping PPE ready where and when it's needed. Constructed from durable, crystal-clear PETG, it offers high visibility, chemical resistance, and a polished, professional appearance that elevates any lab environment.

11 1/4 x 9 3/4 x 3 3/4 inches WHD

Features:

- Dual-sided access provides glove availability from both sides of shared lab areas
- Durable, impact-resistant PETG maintains clarity and strength over time
- Smooth, easy-to-clean surface supports routine disinfection and hygiene compliance



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee