



Element 05 Cart

51504



Designed for care, built for momentum

In a fast-paced clinical setting, the right cart is a critical partner. The Element 05 is built to withstand the day-to-day stresses of its environment as a quiet, secure storage option. Its slim design navigates tight spaces, while durable materials resist stains from disinfectants, chemo, and other harsh chemicals. With smooth casters, locking storage, and a lifetime guarantee, it's a reliable choice for everyday care.

SPACE-OPTIMIZED FOOTPRINT

- Slim width means the cart fits where typical bulky workstations can't, maximizing lab floor utility.

CORROSION-RESISTANT HYBRID CONSTRUCTION

- ACP + HDPE shell resists rust, chipping, and chemical exposure—preserving appearance and function over years.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



[ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



FULL EXTENSION LOCKING DRAWERS

- Three lockable drawers slide fully open, so reagents and tools are effortlessly accessible.



VERSATILE ORGANIZATION

- A combination of side shelves, drawers, and bulk storage, allow for a wide range of supplies to be easily stored.



BULK STORAGE BEHIND CLEAR HINGED DOOR

- Store larger items securely and visually with the lockable bulk storage compartment.



5" LOCKING STABILITY

- Locking 5" casters provide quiet, secure mobility. Handle thresholds and building transitions with ease.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

Element 05 Cart

51504

Drawers & Bulk Storage

- 3 lockable drawers
- 12x6x15 $\frac{3}{8}$ " (WHD)
- Lockable bulk storage area

Acrylic Side Shelves

- 6 Clear Acrylic side shelves
- Hold extra supplies, or dispense gloves/wipes.

Cleanability

- All surfaces are easily cleanable with soap and water, isopropyl alcohol or a hospital-grade cleaner.



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee

ACCESSORIES



PN: 50917
4 Compartment



PN: 50974
6 Compartment



PN: 50975
4 Compartment

