



# Element 07 Cart

## 51610

### MAXIMIZE STORAGE. MINIMIZE FOOTPRINT.

TrippNT's Element 07 Compact Cart delivers serious storage in a design slim enough to fit under counters, tables, and desks. Built with corrosion-resistant materials and smooth-rolling casters, it's the perfect mobile solution for labs where space is limited but efficiency can't be compromised. Four full-extension drawers and a flat workspace top keep supplies organized and accessible, helping your team stay productive and compliant.

### SPACE-SAVING STORAGE

- Four full-extension drawers and a flat worktop deliver maximum organization in a compact footprint.

### EASY MANEUVERABILITY

- Built-in handle and quiet, non-marking casters make transport effortless—even in tight spaces.

### CLEANABLE & DURABLE

- Corrosion-resistant ACP and HDPE construction withstands harsh disinfectants, lab chemicals and daily use.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



#### Built In Handles

- Cleanable one piece top has integrated handle and grooved channel to catch liquids



#### Full-Extension Drawers

- Every inch of storage is visible and accessible. No more wasted time reaching in the dark.



#### 3 inch Casters

- 2 locking and 2 swivel for smooth maneuverability



#### Locking System

- A single key secures all four drawers, protecting sensitive supplies.



# Element 07 Cart

51610

## Compact Design

- Fits under counters, desks, and benches to save valuable floor space.

## Full-Extension Drawers

- Access supplies easily with smooth, soft-close drawers that maximize storage.

## Locking Casters

- Two swivel and two locking 3-inch casters keep the cart steady when needed.

## Durable Construction

- Lightweight aluminum composite and HDPE surfaces resist corrosion and heavy use.

## Flat Worktop

- Creates an instant mobile workstation wherever your team needs it.

## LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Corrosion Resistant



Short Lead Times



Ships Assembled




Made in the U.S.A



100% Lifetime Guarantee

**ASK ABOUT OUR DEMO PROGRAM**

[Try it out at your facility.](#)



## ACCESSORIES



PN: 50917  
4 Compartment



PN: 50974  
6 Compartment



PN: 50975  
4 Compartment



PN: 51998  
Adjustable



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com