



# MRI Core DX Extended Top

## 51817



### Mobilized Storage & Organization for MRI Suites

Streamline your MRI Suite's operations with the MRI Core DX Extended Top Cart—featuring a spacious 29 × 19-inch work surface, secure storage, and built-in glove dispensers. Durable HDPE construction and smooth casters ensure easy movement across shared spaces while meeting strict cleaning and safety protocols

#### EXPANSIVE WORK SURFACE

- Extra-large 29 × 19-inch top gives your team ample space to stage reagents, tools, and documentation—right where they're needed.

#### MR CONDITIONAL CONSTRUCTION

- Built for safety within a 3 Tesla level MRI environment.

#### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



### ASK ABOUT OUR DEMO PROGRAM

[Try it out at your facility.](#)



#### Sliding Door

- Secure items in the bulk storage area with the sliding door.



#### Spacious Drawers

- Three drawers that each measure 11 ½" wide by 15 ¼" deep by 3 ¾" tall.



#### Dispensing Shelf

- Easily access gloves, kimwipes or other supplies.



#### Easy-Grab Handles

- Ergonomic, lightweight, and easy to maneuver.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | [Sales@TrippNT.com](mailto:Sales@TrippNT.com)

# MRI Core DX with Extended Top

51817

## Secure Sliding Door Bulk Storage

- Bulk supplies stay tucked safely behind a sliding door, keeping shared equipment and workflows uncluttered.

## Non-Metal Drawer Slides

- Drawers extend 70% forward with a complete HDPE plastic design, includes a dedent to secure the drawers when closed.



Corrosion Resistant



Short Lead Times




Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



## ACCESSORIES



PN: 50917  
4 Compartment



PN: 50974  
6 Compartment



PN: 50975  
4 Compartment



PN: 51998  
Adjustable