



MRI Core DX Extended Top

51817



BUILT FOR MRI. ENGINEERED FOR EFFICIENCY.

In MRI environments, calm and efficiency are essential. The MRI Core DX with Extended Top supports daily workflow and peace of mind by giving staff a safe, reliable place to store supplies without disruption. It's quiet drawers and non-magnetic construction help maintain a focused, comfortable atmosphere for patients and clinicians alike.

MR CONDITIONAL DESIGN

- Built with minimal metal—only stainless fasteners and aluminum composite drawers.

PREMIUM NON-METAL DRAWER SLIDES

- Durable HDPE slide system with detents offers smooth function without a trace of metal.

EASY ACCESS SIDE WINGS

- Six compartments (two dispensing, four enclosed).
- Ideal for gloves, wipes, and essentials.

EXTENDED WORK SURFACE

- Prep and accomplish more with large extended top.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



[ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



Sliding Door

- Secure items in the bulk storage area with the sliding door.



Non-Metal Slides

- Extend 75% and close smoothly with a durable plastic detent.



Dispensing Shelf

- Easily access gloves, kimwipes or other supplies.



Easy-Grab Handles

- Ergonomic, lightweight, and easy to maneuver.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

MRI Core DX with Extended Top


51817




Premium high density polyethylene surfaces.



High quality metal-free drawer slide with detent!



See Drawer In Action!



- Corrosion Resistant
- Short Lead Times
- Ships Assembled
- Made in the U.S.A
- 100% Lifetime Guarantee

ACCESSORIES

