



# MRI Maui Lab Island

## 51820

### MRI Utility Engineered For Speed And Security

In MRI suites where safety, organization, and speed all must align, the Maui Lab Island Cart delivers. Built of non-ferrous HDPE with a locking tambour door and internal shelf, it gives you secure, visible storage plus a usable workspace. Even fully loaded, it maneuvers easily on 5" twin-wheel casters, letting you stage coils, linens, and consumables close to the magnet. Resistant to chemicals and rigorous cleaning, it maintains a clean, professional look—and with TrippNT's lifetime replacement guarantee, it's a cart that's built to last under pressure.

### 100% MRI-CONDITIONAL MATERIALS

- 99.4% Non-Metal construction. Zero ferrous metals.

### DURABLE, CHEMICAL-RESISTANT SURFACES

- HDPE construction resists stains, disinfectants, and daily wear, preserving aesthetics.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

### [ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



#### LOCKING MECHANISM FOR SECURITY

- Integrated lock ensures spare linens, coils, or sensitive supplies stay in place when unattended.



#### INTERNAL STATIONARY SHELF

- Keeps frequently used items within reach without needing to dig.



#### SLIDING TAMBOR DOOR

- Secures gear with a roll-up door that disappears when open—no swinging panels in tight space.



#### SMOOTH, LARGE CASTER

- 5" twin-wheel casters (2 locking, 2 swivel) provide easy gliding over scan room floors.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | [Sales@TrippNT.com](mailto:Sales@TrippNT.com)

# MRI Maui Lab Island

51820

## SLIM FOOTPRINT WITH AMPLE HEIGHT

- Compact (24.25" × 24.125") footprint that fits into tight corridors—yet tall enough (37.625") for vertical storage.

## KEEP COILS AND PATIENT POSITIONERS CLOSE

- Stay on schedule by creating a flexible and mobile workflows

## LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Corrosion Resistant



Short Lead Times




Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



## RELATED OPTIONS

### MRI Cayman



51821

### MRI Belize



50734