



MRI Cayman Lab Island

51821

High-Capacity Mobile Workbench for Clinical MRI Suites



In a fast-moving MRI suite, where space is tight and every second matters, MRI technologists need a workstation that keeps up with their rhythm. The TrippNT MRI Cayman Lab Island delivers MRI-safe mobility and organized storage right where it's needed most. Built from non-ferrous HDPE and designed for smooth movement in narrow corridors, it provides a stable surface for staging and swapping coils, managing linens, and keeping small consumables within easy reach. With its compact footprint, roll-top access, and easy-clean surfaces, it helps technologists maintain focus, safety, and flow — even on the busiest scan days.

100% MRI-CONDITIONAL MATERIALS


- 99.4% Non-Metal construction. Zero ferrous metals.

DURABLE, CHEMICAL-RESISTANT SURFACES

- HDPE construction resists stains, disinfectants, and daily wear, preserving aesthetics.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



[ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



LOCKING MECHANISM FOR SECURITY

- Integrated lock ensures spare linens, coils, or sensitive supplies stay in place when unattended.



INTERNAL STATIONARY SHELF

- Keeps frequently used items within reach without needing to dig.



SLIDING TAMBOR DOOR

- Secures gear with a roll-up door that disappears when open—no swinging panels in tight space.



SMOOTH, LARGE CASTER

- 5" twin-wheel casters (2 locking, 2 swivel) provide easy gliding over scan room floors.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

HIGH CAPACITY MODULAR STORAGE

- shelf supports coil placement during prep or change-over, keeping equipment stable while freeing both hands for patient setup.

KEEP COILS AND PATIENT POSITIONERS CLOSE

- Stay on schedule by creating a flexible and mobile workflows

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Corrosion Resistant



Short Lead Times







Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



RELATED PRODUCTS



MRI Maui

51820



MRI Belize

50734