



MR Conditional Triple Glove Box Holder





Part No. 51874

Store and access gloves safely within the MRI suite. This triple glove box holder is built from non-metallic, MR Conditional materials—including white PVC and clear acrylic—and is tested to be safe near 3-Tesla magnets. With three compartments and mounting keyholes, it's a dependable, easy-to-clean solution for keeping PPE accessible and compliant inside controlled environments.

11 5/8" x 15 3/8" x 4 5/8" inches WHD

Features:

- **MR Conditional Certified:** Safe for use in environments with static magnetic fields of 3-Tesla or less
- **Triple Box Capacity:** Holds three glove boxes to support high-traffic areas and varied glove sizes
- **Clear Acrylic Front:** Instantly see supply levels and types for fast restocking
- **Wall Mount Ready:** Keyholes on back panel make installation quick and secure



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee