



## BD Sharps Container 5.4 Quart

Part No. 53293

Protect people and maintain compliance with this durable, disposable sharps container, expertly designed to securely hold needles, syringes, lancets, razors, and other sharp objects. Constructed from high-impact, shatter-resistant polypropylene, its punctureproof build ensures reliability in demanding environments. The innovative automatic rotating chamber allows for touchless sharps disposal, reducing accidental contact risks and ensuring hygienic safety. Perfect for warehouses, medical facilities, schools, and rest areas.

12 1/8 x 10 7/8 x 4 3/4 inches WHD



### Features:

- Touchless disposal with automatic rotating chamber
- Puncture-Resistant Design
- Securely contains needles, syringes, lancets, razors, and other sharps
- Fits TrippNT Phlebotomy Pro (51536) and Phlebotomy Pro Mini (53292)

Product Overview



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**