



# Core DX Cart

## Extended Top and Sliding Door

53369



### Everything your lab needs within reach

Streamline your lab's operations with the Core DX Extend Top Cart—featuring a spacious 29 × 19-inch work surface, secure storage, and built-in glove dispensers. Durable HDPE construction and smooth casters ensure easy movement across shared spaces while meeting strict cleaning and safety protocols.

### EXPANSIVE WORK SURFACE

- Extra-large 29 × 19-inch top gives your team ample space to stage reagents, tools, and documentation—right where they're needed.

### Hygienic Durable Construction

- Corrosion-resistant, non-metal HDPE resists harsh disinfectants and laboratory wear—helping you meet compliance without breakdowns.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



### [ASK ABOUT OUR DEMO PROGRAM](#)

Try it out at your facility.



#### Casters

- Ergonomic twin wheel casters, great for traversing bumps and gaps.
- 2 locking, 2 swivel.



#### Spacious Drawers

- Three lockable drawers that each measure 11 1/2" wide by 15 1/4" deep by 3 3/8" tall.



#### Dispensing Shelf

- Easily access gloves, kimwipes or other supplies.



#### Easy-Grab Handles

- Ergonomic, lightweight, and easy to maneuver.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | [Sales@TrippNT.com](mailto:Sales@TrippNT.com)

# Core DX Cart with Extended Top and Sliding Door 53369

## Full Extension Drawer Slides

- Three deep drawers ensure fast, easy access to consumables—helping prevent downtime and streamline inventory control.

## Secure Sliding Door Bulk Storage

- Bulk supplies stay tucked safely behind a sliding door, keeping cold zones, shared equipment, and workflows uncluttered.



Corrosion Resistant



Short Lead Times




Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



## ACCESSORIES



PN: 50917  
4 Compartment



PN: 50974  
6 Compartment



PN: 50975  
4 Compartment



PN: 51998  
Adjustable