



# Core CL Cart Series



## Suture Transport and Storage

The Core CL Cart contains plenty of storage while providing excellent transportation for users on the go. Its corrosion-resistant construction and smooth mobility simplify restocking, cleaning, and inventory control—helping healthcare teams deliver better care and operate more efficiently.

### Suture Storage Shelves

- 6 large built in shelves that can hold up to 36 regular suture boxes at a 45 degree angle to prevent falling.

### Spacious Lockable Drawers

- 3 Drawers measuring 11 ½" x 15 ¼" x 3 ⅜" that can be locked with a key.

### Bulk Storage Area

- A large lockable storage area protected by a clear sliding door.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

### [ASK ABOUT OUR DEMO PROGRAM](#)

Try it out at your facility.



#### Twin-Wheel Casters

- Smooth and easy to move over thresholds and elevator gaps.
- Lockable



#### Spacious Drawers

- Three lockable drawers that each measure 11 ½" wide by 15 ¼" deep by 3 ⅜" tall.



#### Suture Box Shelves

- Holds up to 36 suture boxes for convenient use or restocking



#### Easy-Grab Handles

- Ergonomic, lightweight, and easy to maneuver.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

# Core CL Cart Series

## Twin Casters

- 2 locking and 2 swivel 4" casters

## Easy To Grip Handles

- 3 built-in handles for efficient moving and maneuvering through any space.

## Flat Top Workstation

- Has a flat work surface to write on papers, use a laptop, or to set things on top of.



Corrosion Resistant



Short Lead Times






Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



## COLOR OPTIONS



PN: 53490  
White



PN: 53491  
Almond Beige



PN: 53495  
Global Blue



PN: 53501  
Black



PN: 53592  
Dolphin Gray



PN: 53492  
Silver Metallic

## ACCESSORIES



PN: 50917  
4 Compartment



PN: 50974  
6 Compartment



PN: 50975  
4 Compartment



PN: 51998  
Adjustable