



Core DX Hinged Door & Railtop Series



MOBILE STORAGE BUILT FOR ORGANIZATION, SECURITY, AND DURABILITY

The Core DX Hinged Door Cart with Railtop enhances productivity with smart compartmentalization, secure storage, and fast access to critical supplies. Its corrosion-resistant construction and smooth mobility simplify restocking, cleaning, and inventory control—helping healthcare teams deliver better care and operate more efficiently.

DURABLE & HYGIENIC

- Enjoy lightweight but strong HDPE and ABS construction.
- Easily disinfect and maintain all surfaces.

ORGANIZED STORAGE

- Customize six side shelves for fast access or high-use supplies.
- Maximize storage with three lockable drawers and a bulk storage compartment.

EASY MOBILITY

- Get into tight spaces with smooth casters and built-in easy-grab handles.



ASK ABOUT OUR DEMO PROGRAM

Try it out at your facility.



Twin-Wheel Casters

- Smooth and easy to move over thresholds and elevator gaps.
- Lockable.



Spacious Drawers

- Three lockable drawers that each measure 11 ½" wide by 15 ¼" deep by 3 ¾" tall.



Dispensing Shelf

- Easily access gloves or other supplies.



Railtop Containment

- Prevent supplies, tools, and equipment from falling off the top with the built-in rails.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

Core DX Hinged Door & Railtop

Mobile Supply & Procedure Cart

DRAWERS

- 3 lockable drawers
- Measure 11 1/2" x 15 1/4" x 3 3/8" (WDH)
- Durable ACP interior with HDPE front

BULK STORAGE COMPARTMENT

- Impact-resistant and lockable hinged door
- Measures 12 7/8" x 16" x 15 5/8" (WDH)


CHART & TOOL HOLDER

- 11 3/8" x 1 1/2" x 5" (WDH) main compartment, great for documents or charts
- 1 1/2" x 1 1/2" x 5" (WDH) secondary compartment

CLEANING

- All surfaces can be cleaned with soap and water or hospital-grade cleaner.

-  **Corrosion Resistant**
-  **Short Lead Times**
-  **Ships Assembled**
-  **Made in the U.S.A**
-  **100% Lifetime Guarantee**



COLOR OPTIONS



PN: 53547
White



PN: 53530
Almond Beige



PN: 53534
Global Blue



PN: 53541
Black



PN: 53590
Dolphin Gray



PN: 53531
Silver Metallic

ACCESSORIES



PN: 50917
4 Compartment



PN: 50974
6 Compartment



PN: 50975
4 Compartment



PN: 51998
Adjustable