



# Equipment Carts

## Series



### 14 MODELS. ONE RELIABLE SOLUTION.

From compact 12 inch carts to full-capacity 35 inch models, this family of equipment carts is engineered to mobilize, store, and secure the tools your lab depends on. Every model is built from corrosion-resistant HDPE and PVC, ships fully assembled, and supports up to 400 lbs so you can choose the right size with confidence. Whether your team needs a slim-profile cart for tight spaces or a larger unit with tambour doors and an adjustable shelf, there's a model designed for your exact workflow.

### FULL RANGE OF SIZES

- Fourteen models, from 18 to 30 inches wide and 12 to 35 inches tall, deliver the right fit for every lab environment.

### SECURE STORAGE OPTIONS

- Select models include sliding tambour doors and an adjustable shelf for up high-value equipment.

### HEAVY-DUTY PERFORMANCE

- Each cart supports up to 400 lbs while weighing just 65–85 lbs for easy handling.

### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



#### Mounting Strap & Slots

- Six built-in slots (one strap included) secure equip during movement.



#### Tambour Door

- Prevents Dust and secures supplies.



#### Adjustable Shelf

- Available on taller models

**ASK ABOUT  
OUR DEMO  
PROGRAM**

Try it out at your  
facility.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

# Equipment Carts Series



PN: 53604  
White



PN: 53611  
White



PN: 53623  
White



PN: 53612  
White



PN: 53625  
White



PN: 53614  
White



PN: 53617  
White



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



PN: 53637  
Black



PN: 53640  
Black



PN: 53638  
Black



PN: 53641  
Black



PN: 53639  
Black



PN: 53642  
Black



PN: 53643  
Black



Scan to  
Learn More

Part #	Color	Dimensions (L x W x H) Inches	Tambour Door	Shelf
53604	White	18 x 24 x 12	No	No
53637	Black	18 x 24 x 12	No	No
53611	White	24 x 24 x 12	No	No
53640	Black	24 x 24 x 12	No	No
53623	White	18 x 24 x 12	Yes	Yes
53638	Black	18 x 24 x 24	Yes	Yes
53612	White	24 x 24 x 24	Yes	Yes
53641	Black	24 x 24 x 24	Yes	Yes
53625	White	18 x 24 x 35	Yes	Yes
53639	Black	18 x 24 x 35	Yes	Yes
53614	White	24 x 24 x 35	Yes	Yes
53642	Black	24 x 24 x 35	Yes	Yes
53617	White	30 x 30 x 18	Yes	No
53643	Black	30 x 30 x 18	Yes	No