



Equipment Carts

Series



14 MODELS. ONE RELIABLE SOLUTION.

From compact 12" carts to full-capacity 35" models, this family of equipment carts is engineered to mobilize, store, and secure the tools your lab depends on. Every model is built from corrosion-resistant HDPE and PVC, ships fully assembled, and supports up to 400 lbs so you can choose the right size with confidence. Whether your team needs a slim-profile cart for tight spaces or a larger unit with tambour doors and an adjustable shelf, there's a model designed for your exact workflow.

FULL RANGE OF SIZES

- Fourteen models, from 18" to 30" wide and 12" to 35" tall, deliver the right fit for every lab environment.

SECURE STORAGE OPTIONS

- Select models include sliding tambour doors and an adjustable shelf for up high-value equipment.

HEAVY-DUTY PERFORMANCE

- Each cart supports up to 400 lbs while weighing just 65–85 lbs for easy handling.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Mounting Strap & Slots

- Six built-in slots (one strap included) secure equip during movement.



Tambour Door

- Prevents Dust and secures supplies.



Adjustable Shelf

- Available on taller models

**ASK
ABOUT
OUR DEMO
PROGRAM**

Demo Program

- Try at your facility.



Equipment Carts Series



PN: 53604
White



PN: 53611
White



PN: 53623
White



PN: 53612
White



PN: 53625
White



PN: 53614
White



PN: 53617
White



PN: 53637
Black



PN: 53640
Black



PN: 53638
Black



PN: 53641
Black



PN: 53639
Black



PN: 53642
Black



PN: 53643
Black



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee



Learn more

Part #	Color	Dimensions (LxWxH)	Tambour Door	Shelf
53604	White	18" x 24" x 12"	No	No
53637	Black	18" x 24" x 12"	No	No
53611	White	24" x 24" x 12"	No	No
53640	Black	24" x 24" x 12"	No	No
53623	White	18" x 24" x 12"	Yes	Yes
53638	Black	18" x 24" x 24"	Yes	Yes
53612	White	24" x 24" x 24"	Yes	Yes
53641	Black	24" x 24" x 24"	Yes	Yes
53625	White	18" x 24" x 35"	Yes	Yes
53639	Black	18" x 24" x 35"	Yes	Yes
53614	White	24" x 24" x 35"	Yes	Yes
53642	Black	24" x 24" x 35"	Yes	Yes
53617	White	30" x 30" x 18"	Yes	No
53643	Black	30" x 30" x 18"	Yes	No