



Height Adjustable Double Sided Suture Stand

53944



ORGANIZED SUTURES. FASTER PROCEDURES.

Clinicians need sutures organized, visible, and instantly accessible – yet too often they're buried in drawers or scattered across workstations. This height-adjustable, double-sided suture stand puts up to 48 boxes at eye level, reducing search time, improving posture, and keeping procedures moving smoothly from start to finish.

MAXIMIZED ACCESSIBILITY

- Shelves positioned at a 45° angle let clinicians **see and reach sutures instantly** without digging or bending.

DUAL-SIDE ORGANIZATION

- Work from either side without repositioning the cart. **Double the usable space in the same footprint.**

SMOOTH MOBILITY FOR RAPID DEPLOYMENT

- Casters roll easily over thresholds and between spaces, letting teams **move supplies where they're needed most.**

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



[ASK ABOUT OUR DEMO PROGRAM](#)

[Try it out at your facility.](#)



Height Adjustable

- Raise or lower the cart top quickly and easily with the pneumatic handle



Angled Shelving

- Built to hold suture boxes securely while moving through hallways



Locking Casters

- When the cart is stationary and locked, sutures can be dispensed right from the boxes on the cart



Steering Handle

- The built-in handle allows the cart to move right along with you



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

Height Adjustable Double Sided Suture Stand

53944

ADJUSTABLE HEIGHT CONTROL

- Quick pneumatic lift for easy height changes catering to users of any stature.

SECURE STATIONARY LOCK

- Lock wheels in place for stable suture access during procedures.

EASY PUSH HANDLES

- Located for **ergonomic control** during transport down corridors or OR suites.

OR-READY CAPACITY

- Designed to hold up to **48 standard suture boxes**, keeping your most-used inventory at the point of care.



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee

