



Companion Cart

53981

Adaptable. Mobile. Hospital Grade Plastic.



Every lab professional wants a workspace that adapts as quickly as their workflow. But traditional fixed-height carts limit flexibility and can slow down data entry, analysis, and collaboration. That lack of adaptability can lead to frustration—and even discomfort—during long days at the bench. The Lab Grade Height Adjustable Companion Cart solves this with one-handed pneumatic height control, smooth mobility, and a durable, cleanable HDPE surface. It's the perfect balance of comfort, functionality, and reliability—so your team can stay focused on results, not workspace limitations.

ADJUSTABLE ERGONOMICS

- One-handed pneumatic height control transitions from sitting (28.54") to standing (42.52").
- Encourages better posture and comfort during extended use.

LAB-GRADE CONSTRUCTION

- HDPE top resists stains, chemicals, and moisture for long-term performance.
- Smooth, nonporous surface supports frequent cleaning and disinfection.



Lab Grade Top

- Resistant to chemicals
- Built in Handles
- Velcro strap to secure equipment when moving..



Smooth Mobility

- Open base for easy movement
- 2 locking casters



Ergonomic lever

- Smooth, pneumatic lift control lets users switch from sitting to standing with one hand.



Product Title

Subtitle

MOBILE & STABLE

- Lightweight frame glides easily on four casters; two lock for stability when stationary.
- Ideal for dynamic lab and healthcare workflows.

H-SHAPED BASE DESIGN

- Open base design promotes comfortable foot placement and unobstructed movement.
- Enhances access and ergonomics in tight or shared work areas.

SECURE WORKSURFACE

- Velcro strap keeps laptops or devices secure while the cart is in motion.
- Provides stability during data entry or patient care.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee

