



## 3 Slot Magnet Pipette Holder


Part No. 54001

When your pipettes have a consistent home, your workflow stays smooth. This magnetic holder keeps three pipettes upright and within reach so you can move quickly and stay focused on precision.

6 X 4 1/8 X 5 7/8 inches WHD

### Features:

- Holds any brand of pipettors, 3 slots
- Strong scratch-resistant neodymium magnet mount
- Constructed from polycarbonate
- Stores and organizes Auto pipettes



**Corrosion Resistant**



**Short Lead Times**



**Ships Assembled**



**Made in the U.S.A**



**100% Lifetime Guarantee**