



Vertical Two-Bin Kimwipe Dispenser, Magnet Mount

Part No. 54025


Quick access matters. This dual-bin, magnet-mounted dispenser keeps wipes visible, separate, and within reach — no clutter, no delays. Faster access, cleaner surfaces, less downtime.

12 5/8 x 10 x 4 1/4 inches WHD



Features:

- Holds two boxes of Kimberly Clark number 34133
- Automatically loads box 2 when box 1 is removed
- Made of ultra clear plastic for maximum visibility
- Comes with TrippNT's ultra-strong, soft rubber-coated magnet



Watch the
Video!



Corrosion Resistant



Short Lead Times



Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee