



# MRI Compact Cart

## For the MR Suite & Healthcare



### CONFIDENCE IN MR CARE—COMPACT, CLEANABLE, AND UNMATCHED IN DESIGN.

Purpose-built for MR Conditional environments, the **MRI Compact Cart** gives you peace of mind where safety and space matter most. Its slim footprint, cleanable construction, and secure storage make it a perfect fit under counters or in tight areas of imaging suites. Made from corrosion-resistant plastic with minimal metal, it delivers smart, mobile organization without compromise.

#### COMPACT WITH PURPOSE

- Slim size fits under counters or in tight corners—ideal for maximizing underused space

#### PREMIUM NON-METAL DRAWER SLIDES

- Durable HDPE and TPU slide system offers smooth, soft-close function without a trace of metal.

#### MR CONDITIONAL CONSTRUCTION

- Safe for MR suites with minimal metal—only stainless fasteners and aluminum composite drawer fronts.

#### CLEANABLE & CORROSION-RESISTANT

- Smooth HDPE and ACP surfaces are clinically ready, resisting damage and withstanding harsh disinfectants.

#### LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

### ASK ABOUT OUR DEMO PROGRAM

Try it out at your facility.



#### Built-In Top Containment

- Use the top for secure & fast access.



#### Lockable Casters

- 2" lockable casters let you move and secure the MRI Compact Cart with ease.



#### Premium Soft-Close Drawers

- Enjoy our all-plastic design without compromising quality!



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

## MRI Compact Cart


For the MR Suite & Healthcare

**With 99.4% Non-Ferrous Construction,  
No Other MR Cart Comes Close.**

- ✓ No metal-body construction,
- ✓ No metal drawer slides,
- ✓ And no loss of functionality!



Premium high density  
polyethylene surfaces.



Soft-close mechanism in  
every drawer slide!

## ACCESSORIES



PN: 50917  
4 Compartment



PN: 50974  
6 Compartment



PN: 50975  
4 Compartment



PN: 51998  
Adjustable