



LabFlow Hood Cart

Underhood Tray Organization



ORGANIZED STORAGE RIGHT UNDER YOUR HOOD

The **LabFlow Hood Cart** keeps essential supplies directly under your biosafety cabinet or fume hood, organized, accessible, and ready to go. Designed for high-efficiency lab workflows, it eliminates clutter and frustration by giving every tray a designated place. Built entirely from cleanable HDPE and featuring smooth, locking casters, it rolls in easily and stays put when you need it to.

POINT-OF-USE HOOD STORAGE

- At under 27" tall, keep your supplies and tools accessible directly beneath the work surface.
- Ideal for improving efficiency and reducing workspace mess.

MODULAR TRAY SYSTEM

- Holds standard LabFlow trays in multiple depths for flexible storage.
- Trays transfer easily between other LabFlow products, including cold storage.

HIGH VISIBILITY ACCESS

- Open-front design makes tray contents visible at a glance.
- Speeds up supply retrieval and supports fast-paced workflows.

LAB-GRADE DURABILITY

- HDPE enclosure resists damage, staining, and harsh chemicals.
- Aluminum rails resist corrosion and are easy to clean.
- Wipes down quickly to meet lab hygiene standards.

LIFETIME GUARANTEE

- Work stress-free, knowing your cart is covered for life.

ASK ABOUT OUR
DEMO
PROGRAM

Try it out at your facility.



Built-In Top Containment

- Use the top for secure & fast tray access.



Fits LabFlow Trays

- Receives LabFlow Trays of any depth, size, or material.



Flexible Rail System

- Six levels of rails let you accommodate LabFlow Trays of various sizes.



Low Profile Casters

- 2" lockable casters let you roll the LabFlow Hood Cart in or out with ease.



LabFlow Hood Cart

Underhood Tray Organization



MULTIPLE TRAY OPTIONS



Learn about
LabFlow



Material	Cart with 6" Trays	Cart with 4" Trays	Cart with 2.5" Trays
Polycarbonate	PN 54096	PN 54097	PN 54098
Polypropylene	PN 54099	PN 54100	PN 54101
No Trays	54089		



- Corrosion Resistant
- Short Lead Times
- Ships Assembled
- Made in the U.S.A
- 100% Lifetime Guarantee

ADDITIONAL ACCESSORY TRAYS


Polycarbonate

Polypropylene



PN	Description
101172	Polycarbonate 2.5"
101167	Polycarbonate 4"
101120	Polycarbonate 6"

PN	Description
101171	Polypropylene 2.5"
101126	Polypropylene 4"
101170	Polypropylene 6"



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com