



LabFlow Storage Rack

For Laboratories and Cold Storage



CONQUER YOUR STORAGE CLUTTER

The **LabFlow Storage Rack** doubles cold storage capacity by reclaiming vertical space with 25 rail levels and secure trays. It prevents spills, protects samples, and keeps inventory visible for fast audits. Clear, color-coded trays reduce searching and cross-contamination, while durable aluminum and locking casters ensure long-term reliability and mobility.

COMPATIBLE WITH COLD STORAGE

- Performs in cold environments like 4°C walk-ins without warping or corroding.
- Keeps supplies off open shelves and sorted for faster retrieval.

MODULAR TRAY SYSTEM

- 25 rail levels use vertical space more efficiently than wire shelving.
- Holds any LabFlow standard tray.

VISIBILITY AT A GLANCE

- Clear trays and open-side rails let you see contents instantly.
- Makes inventory checks faster and improves accountability in shared storage areas.

LAB-GRADE DURABILITY

- High quality welded and riveted aluminum construction.
- Will hold up to even rough treatment and will never corrode.

LIFETIME GUARANTEE

- Work stress-free, knowing your storage rack is covered for life.



Corrosion Resistant



Short Lead Times




Ships Assembled



Made in the U.S.A



100% Lifetime Guarantee




Watch the
Product Overview
Video!



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com

LabFlow Storage Rack

For Laboratories and Cold Storage



Manual Gate Lock

- The front and back each have a built-in manual gate lock for securing trays inside the rack during transportation.

Material	Cart with 6" Trays	Cart with 4" Trays	Cart with NO Trays
Polycarbonate	PN 54110	PN 54111	PN 54092
Polypropylene	PN 54108	PN 54109	



Lockable Casters

- The two front casters are lockable for securing the rack in place.

Dimensions: 17 ³/₄" x 69 ³/₈" x 21 ⁵/₈" (WHD)

COMPATIBLE WITH LABFLOW TRAYS

Polycarbonate



Polypropylene



PN	Description	Max Qty
101172	Polycarbonate 2.5"	13
101167	Polycarbonate 4"	13
101120	Polycarbonate 6"	9

PN	Description	Max Qty
101171	Polypropylene 2.5"	13
101126	Polypropylene 4"	13
101170	Polypropylene 6"	9

**ASK ABOUT OUR
DEMO
PROGRAM**

Try it out at your facility.



10991 N Airworld Drive, Kansas City, MO, USA | 1 (800) 874-7768 | Sales@TrippNT.com