



# Single Hook for Adjustable Surgical Helmet Rack





Part No. 53942

TrippNT's single 2 prong hook can be purchased to replace or add additional hooks to your 24 inch or 48 inch adjustable surgical helmet rack. Constructed from non-metal white high-density polyethylene (HDPE), this product is corrosion resistant and easily sterilized. Each hook has a 50 pound capacity.



## Features:

- *Single hook with 2 prongs*
- *Hand tightened thumb screws*
- *Made from high-density polyethylene*
- *Replaces or adds additional hooks to rack*
- *50 pound capacity per hook*



Rev: 53942.1



**Ships Next Day**



**Made in the U.S.A**



**Corrosion Resistant**



**100% Lifetime Guarantee**